

P:\SW-Ftek\Timo en Cleo\Timo 6034

Signed	KA	Date	02-03-2016
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Overview parts Timo


Copyright reserved by law
Unit of measurement: mm

Lopital Nederland bv

Type nr: 54006034

Timo 54006034

Part nr.	QTY.	Art.nr.	Part
1	1	70203501	Frame
2	1	70203509	Backrest
3	1	69608024	Footrest
4	4	70519928	Castor Ø 125 mm braked
5	1	70602007	Seat
6	1	70565001	Armrest
7	4	70602006	Snaploc
8	2		Cylinder head screw M8x12
9	4		Cap castor
10	4		Tube cap
11	4		Countersunk screw M6x10
12	4		Capring M4
13	4		Chipboard screw M4x30
14	4		Screw cap M4



P:\SW-Ftek\Timo en Cleo\Timo 6034

Signed	KA	Date	02-03-2016
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Overview parts Timo

Copyright reserved
by law

Unit of measure-
ment: mm

Lopital Nederland bv

Type nr: 54006034