



P:\SW-Ftek\Timo en Cleo\Timo 6021\

Signed	KA	Date	01-03-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

Copyright reserved
by law

Unit of measure-
ment: mm

Overview parts Timo

Lopital Nederland bv

Type nr: 54006021

Timo 54006021

Part nr.	QTY.	Art.nr.	Part
1	1	70602101	Frame
2	1	70203509	Backrest
3	1	70602105	Seat
4	1	70565001	Armrest
5	4	70602006	Snaploc
6	1	70510050	Guide for bedpan set
7	1	69501104	Bedpan + lid
8	4		Insert Ø25
9	4		Tube cap Ø25
10	2	75005640	Cylinder head screw M8x16 Super Duplex 4529
11	4		Countersunk screw M6x10
12	4		Flat washer M5
13	4		Cylinder head screw M5x8
14	4		Capring M4
15	4		Chipboard screw 4,2x30
16	4		Screw cap M4



P:\SW-Tek\Timo en Cleo\Timo 6021\

Signed	KA	Date	01-03-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

Overview parts Timo

Copyright reserved
by law

Unit of measure-
ment: mm

Lopital Nederland bv

Type nr: 54006021