



P:\SW-Tek\Timo en Cleo\Cleo 6050\

Signed	KA	Date	29-02-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

Overview parts Cleo

Copyright reserved by law
Unit of measurement: mm

Lopital Nederland bv

Type nr: 54006050

Cleo 54006050

Part nr.	QTY.	Art.nr.	Part
1	1	70605001	Frame
2	1	70203509	Backrest
3	2	70300504	Wheel
4	1	70300505	Brake lever (set)
5	2	70300506	Castor Ø125 mm unbraked
6	1	69608024	Footrest
7	1	70565001	Armrest
8	1	70602105	Seat
9	4	70602006	Snaploc
10	1	70510050	Guide for bedpan set
11	1	69501104	Bedpan + lid
12	2		Cap castor
13	6		Insert Ø25
14	2		Hex bolt M12x80
15	2		Washer M12
16	2	75005640	Cylinder head screw M8x16 Super Duplex 4529
17	4		Countersunk screw M6x10
18	4		Flat washer M5
19	4		Cylinder head screw M5x8
20	4		Capring M4
21	4		Chipboard screw 4,2x30
22	4		Screw cap M4



Signed	KA	Date	29-02-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

Copyright reserved
by law

Unit of measure-
ment: mm

Overview parts Cleo

Lopital Nederland bv

Type nr: 54006050