

Z:\SW-Tek\Luxal 2200\


Signed	PS	Date	03-10-2007
Revised	4 /ES	Date	25-02-2013
Revised	5 /KA	Date	24-02-2015

# Overview parts Luxal

Copyright reserved by law	Unit for measurement: mm
---------------------------	--------------------------

Lopital Nederland bv

Type nr: 61002200



Z:\SW-Tek\Luxal 2200\

Signed	PS	Date	03-10-2007
Revised	4 /ES	Date	25-02-2013
Revised	5 /KA	Date	24-02-2015

# Overview parts Luxal Base frame

Copyrights reserved by law

Unit for measurement: mm

Lopital Nederland bv

Type nr: 61002200

## Luxal 61002200

Part nr.	QTY.	Art.nr.	Part
1	1	70220001	Chassis
2	2	70220004	Side rail
3	1	70230007	Frame
4	4	70220026	Wheel bolt
5	1	70220020	Axis long
6	1	70220030	Axis short
7	1	70220021	Coupling rod
8	1	70220022	Coupling rod
9	1	70220019	Hinge
10	1	70210003	Adjusting nut
11	1	70220025	Battery plate
12	1	70220018	Hinge split
13	2	70220017	Axis
14	2	70220029	Bronze bush
15	1	70210512	Mountingplate wall charger
16	1	70210516	Charger
17	2	70220006	Support wheelcover
18	4	70220007	Hook wheelcover
19	1	70220028	Hose clip
20	1	70220016	Mattress
21	1	70210048	Pillow
22	4	70220040	Plug for buffer role
23	2	70220052	Theaded strip clamp
24	4	70220051	Theaded strip Bolt clamp
25	8	70220023	Clamping plate
26	1	70220035	Wheel cover closed
27	1	70220015	Wheel cover open
28	2	70210039	Lock
29	2	70210041	Leaf spring incl. cover
30	1	70210056	Clamping ring
31	1	70210055	Outlet
32	1	70210043	Clamping disk
33	1	70220014	Web plate
34	1	70220033	Brake pedal right or left
35	2	70220009	Mounting bracket
36	2	70220008	Actuator BL4
37	4	70220041	Buffer roll
38	1	70210058	Pakking ring
39	1	70210054	Hose (complete)
40	1	70210511	Mounting plate
41	1	70210518	Batterypack
42	1	70220012	Control box

Signed	PS	Date	03-10-2007
Revised	4 /ES	Date	25-02-2013
Revised	5 /KA	Date	24-02-2015

## Overview parts Luxal

Copyrights reserved  
by law


Unit for measure-  
ment: mm

Lopital Nederland bv

Type nr: 61002200

## Luxal 61002200

Part nr.	QTY.	Art.nr.	Part
43	1	70220013	Hand control
44	1	28021090012	Emergency guard
45	1	70220010	Cable
46	1	70220011	Cable
47	6	80220065	Strap with buckle
48	4	70220005	Castor set
49	24	70220027	Threaded bush
50	2	70220005	Carrier
51	4		Tube cap
52	2		Slide bearing
53	4		Washer
54	2		Slide bearing
55	2		Cylinder head screw
56	1		Hex bolt
57	1		Hex bolt
58	2		Slide bearing
59	4		Retaining ring
60	4		Washer
61	2		Retaining ring
62	4		Slide bearing
63	1		Round head screw
64	4		Hex bolt
65	4		Tube cap
66	20		Round head screw
67	4		Washer
68	4		Rivet
69	46		Bolt
70	11		Washer
71	12		Spring washer
72	8		Lock nut
73	6		Hex bolt
74	8		Hex bolt
75	12		Locknut
76	12		Washer
77	5		Locknut



Z:\SW-Tek\Luxal 2200\

Signed	PS	Date	03-10-2007
Revised	4 /ES	Date	25-02-2013
Revised	5 /KA	Date	24-02-2015

## Overview parts Luxal

Copyrights reserved  
by law

Unit for measure-  
ment: mm

Lopital Nederland bv

Type nr: 61002200