



Z:\SW-Tek\Peilkaan lijn\4050\

Signed	KA	Date	22-01-2015
Revised	- /	Date	-
Revised	- /	Date	-

Copyrights reserved by law	Unit measurement: mm
----------------------------	----------------------

# Overview parts Amfora

Lopital Nederland bv

Type nr: 6200 4050/4055

## Amfora 62004050 / 62004055

Part nr.	QTY.	Art.nr.	Part
1	1	70405001 70405501	Frame L=160 Frame L=190
2	2	67904120 67904125	Side rail complete L=160 Side rail complete L=190
3	1	67908070 67908068	Drip tray L=160 + Hose Drip tray L=190 + Hose
4	4	70519928	Castor Ø 125 mm braked
5	1	70410003 70410501	Inlay frame backrest L=160 Inlay frame backrest L=190
6	1	70410001	Adjusting backrest
7	1	70410004	Clamp high
8	2	70401012 70401512	Lying surface 160 Lying surface 190
9	2	70410002	Tube backrest adjusting
10	1	70404001	Clamp low
11	4		Cap castor
12	10		Washer
13	4		Lock nut
14	4		Nut cap
15	2		Wingbolt
16	2		Tube cap



Z:\SW-Tek\Peilkaan lijn\4050\

Signed	KA	Date	22-01-2015
Revised	- /	Date	-
Revised	- /	Date	-

Copyrights reserved by law	Unit measurement: mm
----------------------------	----------------------

## Overview parts Amfora

Lopital Nederland bv

Type nr: 6200 4050/4055