



P:\SW-Ftek\Timo en Cleo\Cleo 6040\

Signed	KA	Date	29-02-2016
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Overview parts Cleo

Copyright reserved by law

Unit of measurement: mm

Lopital Nederland bv

Type nr: 54006040

Cleo 54006040

Part nr.	QTY.	Art.nr.	Part
1	1	70605001	Frame
2	1	70203509	Backrest
3	2	70300504	Wheel
4	1	70300505	Brake lever (set)
5	2	70300506	Castor Ø 125 mm
6	2	69608024	Footrest
7	1	70565001	Armrest
8	1	70602007	Seat
9	4	70602006	Snaploc
10	2		Cap castor
11	6		Tube cap
12	2		Hex bolt M12x80
13	2		Washer M12
14	2		Cylinder head screw M8x12
15	4		Countersunk screw M6x10
16	4		Chipboard screw M4x30
17	4		Capring M4
18	4		Screw cap M4



P:\SW-Ftek\Timo en Cleo\Cleo 6040\

Signed	KA	Date	29-02-2016
Revised	- /KA	Date	-
Revised	- /KA	Date	-

Overview parts Cleo

Copyright reserved
by law

Unit of measure-
ment: mm

Lopital Nederland bv

Type nr: 54006040